

EE122B

Introduction to Biomedical Electronics

Tissue Impedance

Greg Kovacs, Laurent Giovangrandi

Department of Electrical Engineering
Stanford University

Content

Review of tissue impedance concepts and the use of impedance measurement in the determination of physiological parameters.

- Tissue Impedance Characteristics
- Impedance Measurement Techniques
- Applications
 - Impedance plethysmography (Respiration)
 - Impedance cardiography (ICG)
 - Tissue characterization (Body Composition Analyzers)
 - Electrical Impedance Tomography (EIT)

Science

Tissue Impedance Characteristics

A Few Basics First

- Tissue bioimpedance includes electrical conduction from both a conduction and dielectric point of view, and focuses on two main parameters:
 - Conductivity (charge transport)
 - Permittivity (charge storage)

- Current transport in a lossy medium is defined by:


- σ : conductivity
- ε : permittivity

$$J = J_c + J_d = \sigma E - i\omega\varepsilon E$$

Or: $J = -i\omega\hat{\varepsilon}E$, with $\hat{\varepsilon} = \varepsilon + i\frac{\sigma}{\omega}$ (complex permittivity)


- Permittivity is usually frequency-dependent, due to *relaxation* effects (dipoles can follow field change only so fast), leading to so-called *dispersions* (frequency response peaks) in the frequency domain.

Tissue Impedance Model




At low frequencies (DC-kHz), the current goes around the cells. Higher, it flows through the cells.

 HF
 LF



Simplified equivalent model for conductivity



Dielectric dispersion regimes:

α (mHz-kHz): counter-ion effect around membrane glyocalix, ionic diffusion, intracellular structures.

β (kHz-100MHz): Maxwell-Wagner effect, membrane capacitance, protein response.

γ (GHz): polarization of water molecules

Source: Schwan, Annals of Biomedical Engineering, 16, 1988

Tissue Conductivity Examples




Figure 4.10 Dispersion of muscle tissue. (a) transversal; (b) longitudinal. Source: From Gabriel et al. (1996b) by permission.




Figure 4.11 Dispersion of brain tissue. Source: From Gabriel et al. (1996b) by permission.




Figure 4.12 Dispersion of adipose tissue. Source: From Gabriel et al. (1996b) by permission.


Tissue Conductivity Examples

TABLE 4.1 Tissue Conductivity, LF and 1 MHz

Tissue	σ (S/m) 1 Hz–10 kHz	σ (S/m) ca 1 MHz	$\phi_{\max}$ @<10 MHz	anisotropy
Human skin, dry	10^{-7}	10^{-4}	80°	?
Human skin, wet	10^{-5}	10^{-4}	30°	?
Bone	0.005–0.06		20°	Strong
Fat	0.02–0.05	0.02–0.05	3°	Small
Lung	0.05–0.4	0.1–0.6	15°	Local
Brain (grey matter)	0.03–0.4	0.15	15°	Small
Brain (white matter)	0.03–0.3			Strong
Liver	0.2	0.3	5°	?
Muscle	0.05–0.4	0.6	30°	Strong
Whole blood	0.7	0.7	20°	Flow dependent
Urine	0.5–2.6	0.5–2.6	0°	0
CSF	1.6	1.6	0°	0
Saline, 0.9%, 20°C	1.3	1.3	0°	0
Saline, 0.9%, 37°C	2	2	0°	0
Seawater	5	5	0°	0

A key research question is whether or not malignancies can be detected via impedance. (e.g., <http://www.ncbi.nlm.nih.gov/pubmed/8058021>)

Tissue Resistivity in Thorax



Tissue	Resistivity [Ohm·m]
Blood	1.6
Heart muscle	2.5 parallel 5.6 normal
Skeletal muscle	1.9 parallel 13.2 normal
Lungs	20
Fat	25
Bone	177


Body Impedance Examples

At low frequencies (100Hz-100kHz):

- Whole-body (hand-foot): 500 Ω
- Limb: 100-200 Ω
- Trunk: 50-100 Ω
 - Variation due to respiration:
1-10 Ω
 - Variation due to blood flow:
100-500 m Ω



Requires m Ω resolution !!
(Kelvin connections, typically)



Technology

Impedance Measurements

Impedance Measurement Methods

- Bridges (LCR meter)

Adjustment of variable element to match unknown impedance.

- I - V method

Direct measurement of current/voltage characteristic.


– Lock-in amplification

- Time domain spectroscopy


Broadband excitation recorded in time-domain, analyzed in frequency domain.

- Time domain transmissometry

Model-based current-voltage fitting.




Is It Safe?

$$Z = \frac{V}{I}$$




Of course, isolation is necessary.

- Impedance measurement implies injection of a current.
Is it safe?
- Yes... but must be careful.
 - Keep excitation above 1 kHz!
Perception threshold increases with frequency.
 - Keep current below 4 mA RMS
Generally accepted safe limit.





Electrode Configurations

Preferably, a four-electrode configuration (Kelvin) is used to measure tissue impedance, in order to reduce contributions of the electrodes ($1\text{-}10\text{ k}\Omega$) on the measured impedance (usually $10\text{-}1000\text{ }\Omega$):




$$Z_t = \frac{V}{I} - 2Z_e \quad \text{Unknown}$$



Direct: $Z_t = \frac{V}{I}$


Lock-in Amplification

Lock-in amplifiers act as extreme narrow-band filters which selectively amplify a signal of a given frequency.



Lock-in Amplification

By using two reference signals with a 90° phase-shift (sine/cosine), the lock-in amplifier can separate amplitude and phase:



$$R = \sqrt{X^2 + Y^2} = V_{sig} \quad \varphi = \arctan\left(\frac{X}{Y}\right)$$

Squarewave Lock-in Amplification



- Can do this also using any waveform (e.g., squarewave), whereby each harmonic is also heterodyned (more correctly, homodyned) to DC.
- The, of course, one is “sampling” a bigger fraction of the impedance of the target system and lumping the results together.
- For basic research, single-frequency lock-in approaches are generally used.
- For low-cost consumer applications, squarewave techniques are often used because the “multiplier” is simply a switch driven digitally to pass or block the signal.

Phase-Sensitive Demodulator

- The heart of the lock-in amplifier is the synchronous rectifier, consisting of a phase-sensitive detector (PSD) and a low pass-filter.
- Phase-sensitive detector:
 - Multiplier (analog or digital)
 - Switched gain stage (+1/-1 gain)
 - Diode bridge (double balanced mixer)


Typically used for RF applications

Phase-Sensitive Demodulator



Lock-in Amplifier

Synchronous demodulator based on switched-cap block (Linear Technologies LT1043)




Note the zero-crossing detector to compensate for potential phase shift


Analog Devices AD630 Balanced Modulator/Demodulator

Source: Analog Devices, Application Note AN-683

AD630 Detail



AD633: True Multiplier



00786-023

Commercial Lock-In Amplifiers

Stanford Research Systems (SRS) SR-830

<http://www.thinksrs.com/index.htm>

Digital multiplication!



Zurich Instruments HF2LI Lock-in amplifier

<http://www.zhinst.com/>




Applications

Respiration Measurements


Impedance Plethysmography

- Impedance plethysmography (“measurement of volume”) can be readily applied to respiration monitoring.
- Intake of air in lungs increases thoracic resistance.
- 1-10 Ω modulation, 0.01-0.5 Hz bandwidth.
- Can be calibrated to measure tidal volume.
- Sensitive to motion artifacts.





Ambulatory Systems

NASA/Stanford LifeGuard CPOD
Multiparameter physiologic monitor



getemed VitaGuard® VG2100
Heart and Apnea monitor



Applications


Impedance Cardiography

Impedance Cardiography


- Impedance cardiography (ICG) aims at evaluating hemodynamic parameters:
 - Heart rate
 - Stroke volume (amount of blood ejected at each heart beat), cardiac output
 - Timings (Left Ventricular Ejection Time, Pre-ejection Period,...)
 - Fluid status (Z_0).
- 10-200 m Ω modulation, DC-30 Hz bandwidth.
- Typically, 50-100 kHz/1-4 mA RMS current.



ICG Electrode Setup




Original band-type electrodes



Current discrete gel-type electrodes (8 total)

Source: Bioelectromagnetism, Malmivuo & Plonsey, <http://www.bem.fi> & <http://www.cardiodynamics.com/>

ICG Waveform



Points:

- A:** Contraction of atrium
- B:** Opening of aortic valve
- C:** Maximum flow from left ventricle into aorta
- X:** Closing of aortic valve
- Y:** Closing of pulmonary valve
- O:** Widest open of mitral valve

Signals:

- Z:** Impedance
- dZ/dt :** derivative of impedance
- ECG:** Electrocardiogram
- PCG:** Phonocardiogram

Cardiac Output Derivation


- Challenge: deriving the stroke volume (SV) & cardiac output (CO) from changes in thoracic impedance.

- Nyboer: $\Delta v = \Delta Z \rho \left(\frac{L}{Z_0} \right)^2 \quad [\text{m}^3]$
(1950)

- Kubicek: $SV = \rho \left(\frac{L}{Z_0} \right)^2 \left| \frac{dZ}{dt} \right|_{\min} \cdot t_e \quad [\text{m}^3]$
(1966)
 $CO = SV \cdot HR$

SV stroke volume [ml]
 ρ blood resistivity [$\Omega \cdot \text{cm}$], typically 160 $\Omega \cdot \text{cm}$
 L distance between sensing electrodes
 Z_0 mean impedance of the thorax [Ω]
 $\left| \frac{dZ}{dt} \right|_{\min}$ absolute value of the maximum deviation of the first derivative signal during systole [Ω/s]
 t_e ejection time (=LVET) [s]
 CO cardiac output [l/min]
 HR heart rate [bpm]

Many empirical refinements have been made, and are still being made...



Commercial Systems

CardioDynamics/SonoSite



FDA-approved

Niccomo™/Medis
(Now CardioDynamics/SonoSite)



CardioScreen/Medis




Applications

Body Composition Analysis

Body Composition Analysis

- Estimation of body composition based on whole-body or segmental impedance:
 - Total Body Water (TBW)
 - Extracellular Water (ECW)/Intracellular Water (ICW) balance
 - Fat-Free Mass (FFM, or Body Cell Mass or Muscle Mass)
 - Fat Mass
- Various Configurations:
 - Single frequency (typ. 50 kHz)
 - Multiple frequencies (5-500 kHz)
 - Whole body
 - Segmental




Body Composition Analysis

- Whole-body, single frequency (SF-BIA – most common)
 - TBW from resistance, FFM & FM derived:

$$TBW = k_1 \frac{H^2}{R} + k_2 W + k_3$$

- Phase sometimes used as index of nutrition.
 - Typical errors of 3-4 l for TBW and 3-4 kg for FFM).
 - Ok for home, but not for clinical use!
- Refinements (multiple frequencies, segments) help selectively, but the method remains very sensitive to *calibration* and *fluid abnormalities*.




Omron scale



BioScan 916, Maltron
<http://www.maltronint.com>


Body Composition Example

- Wiithings Scale data, fat mass vs time, ≈ 150 lb male.



Bioelectrical Impedance Vector Analysis

- 'Black-box' approach looking only at patterns in height-normalized resistance/reactance vector.
- No regression, no assumption about models.
- Useful for trending.





After: Kyle, Clinical Nutrition, 23, 2004

Applications


Electrical Impedance Tomography

Electrical Impedance Tomography

- By using many electrodes to take independent measurements, tissue impedance can be imaged.
- Typically 8-16 frequencies used, spanning 10kHz-1MHz.
- 16-128 electrodes, currents in the milliampere range.
- Resolution currently limited ($\approx$ centimeters).
- High sensitivity to anisotropy of tissue in these frequencies.
- Limited clinical use so far.




EIT: Reconstruction Challenges




Phantom Chest

Problem is ill-defined: many degree of freedom, few inputs (inverse conductivity problem)




'D-bar' algorithm




'NOSER' algorithm

EIT Examples

Aide to breast cancer detection




T-Scan 2000 from Siemens (no longer available)
Source: Malich *et al.*, Eur. Radiology, 2000



ITM, LD Technologies (256 electrodes)
<http://www.ldteck.com/>

Chest imaging





Kerrouche, *et al.*, Physiological Measurement, 21, 2001

Applications

Other Applications of Impedance Measurements

Galvanic Skin Response (GSR)



- Changes in skin resistance associated with the modulation of sweat (eccrine) glands activity, also known as *Electrodermal Response* (EDR).
- Measured as *Galvanic Skin Resistance* (GSR) or *Conductance* (GSC).
- Most pronounced on plantar/palmar (glabrous) regions.
- Measured at DC (electrode polarization!) or low-frequency AC.
- Indicator of the psychological state (relaxed, nervous, anxious).
- Used in biofeedback, polygraph.



Sources: Webster & <http://www.iworx.com>

Blood Flow in Limbs

- Venous Occlusion Plethysmography
 - Measure arterial and venous flows
 - e.g., diagnosis of deep vein thrombosis
- Arterial Pulse Wave
 - Arteriosclerosis, stenosis,...



Source: <http://www.medis-de.com>

Lab this week

Respiration and Impedance Cardiography (ICG)